


CDMAM Analyser

Accurate and easy quality assurance

Introduction

- In combination with Artinis CDMAM 3.4 Phantom
- Easy to use
- A necessity for maintaining quality
- Generates quality reports




Characteristics


- Automatic scoring diagram
- Contrast Detail curve
- Reference Contrast Detail curve
- Image quality figure
- Based on EUREF CDCOM program
- User definable number of images

Formats

- DICOM Format
- Runs on Windows XP® platforms
- Free USB port needed



Contrast Detail Curve: $KQFm = 102.57$, Number of analysed images: 8, Detection rate(%): 62.5
Reference Curve: $KQFm = 112.40$, Number of analysed images: 8, Detection rate(%): 62.5



Company Information

ARTINIS MEDICAL SYSTEMS BV
St. Walburg 4
6671 AS Zetten
The Netherlands

Tel: +31 488 422 710
Fax: +31 488 422 829

Email: askforinfo@artinis.com
Web: www.artinis.com

