



TO PAN

DENTAL

TO PAN

A test object that allows the complex focal trough of the panoramic x-ray unit to be quantified. This allows full assessment of parameters affecting image quality. Can be used in film and digital systems.

Allows:

- Quantification of the focal trough and image layer
- Accurate determination of the bite offset
- Quantification of the symmetry of the scan
- Measurement of the incident beams angles
- Simple visualisation of faults



Radiograph of TO PAN

For more information or to find your regional distributor please visit our website

www.leadstestobjects.com | tel. +44(0)1423 321102

Copyright © 2011 Leeds Test Objects Ltd | SpecSh 307a 16/02/2012