

MiRo House Becklands Close Boroughbridge North Yorkshire YO51 9NR UNITED KINGDOM

tel: +44(0)1423 320 007



MagIQ Oncology MRI Phantom

MagIQ Oncology is designed for image quality assessment in the diagnosis, treatment planning and therapy of cancer (especially breast and prostate tumours). Comprises cylindrical phantom which nests into a body section.

Measure:

- Homogeneity
- Contrast Resolution
- Slice Thickness
- MTF / Line Pairs
- · Spherical Gland Volume
- · Geometric Distortion
- Silicon Suppression Techniques
- Fat Saturation Techniques
- Detachable Section for T1 & T2 values and DCE
- Detachable Section for Gddoped Water

External dimensions:

Cylindrical phantom 142mm diameter 180mm deep

Body section 180 x 200 x 300mm






MiRo House Becklands Close Boroughbridge North Yorkshire YO51 9NR UNITED KINGDOM

tel: +44(0)1423 320 007



MagIQ Oncology MRI Phantom







MiRo House Becklands Close Boroughbridge North Yorkshire YO51 9NR UNITED KINGDOM

tel: +44(0)1423 320 007



MagIQ Oncology MRI Phantom



