



RT EPID

Radiotherapy Phantom

RT EPID Features

- Planar MV Imaging QA
- Support for all EPID devices, including:
 - Varian aS500 / aS1000 / aS1200
 - Elekta iView
 - Siemens EPID
- All TG-142 recommended MV imaging tests
- Simple, flexible setup
- Innovative design
- Low cost solution

Routine EPID QA

- Geometric accuracy
- EPID signal isotropy
- Resolution
- Low contrast detectability
- Uniformity and noise
- Automated tracking, trending, and reporting

General Features

- 25cm x 25 cm (tests large portion of panel)
- Lightweight
- Alignment marks for room lasers
- Place directly on EPID surface
- Calibrated for any EPID imager



www.leedstestobjects.com